

Lab Equipments are also evolving like your smart phone in hand



Now, it's

“SMART is the right way to do something in the Lab as well.”

Some tens of equipment are required in the lab for any kind of chemical or biological experiment.


A fact is that a scientist has to walk around in the lab periodically and monitor every equipment to watch if the experimental procedure is going well without any problem.

In spite of the limited manpower, time and space, scientists have a common mission to get the best result in their experiments, and it's not a simple issue to incorporate a certain innovative way to improve laboratory productivity.

Now, DAIHAN Scientific is going to take an initiative to change existing situation in the labs. So called, **Smart-Lab™ Equipment** is being introduced now here for the first time in the world, which will definitely ***change the way to work in the labs.***

and smart TV at home

the time to innovate your lab with Smart-Lab™ Equipments



“Remote Laboratory” - No limitation on Time and Place

- It's possible to monitor and control Wisdom Smart-Lab™ Equipments through App or Web at any time and anywhere you are.

“Intuitive Use” - Easiest way to monitor and control


- You would have Innovative User Experience on the full touch-screen LCD controller
- All the equipments in a lab, or, even in multiple labs, can be controlled at your fingertips using Smart phone or Laptop.

“Secure Laboratory” - Self-diagnostics & Real-time notification

- Based on Self-Diagnostic System, any operational error of the equipment is noticed to the user and to the Service Center, automatically.

Maximum Convenience by Advanced Digital Technology

Smart-Lab™ Equipments



Like your Smart Phone in hand, and Smart TV at Home,
Now, the Various Kinds of Equipment in Your Lab can be SMART!!

Make your Works and Laboratories SMART
with DAIHAN Smart-Lab™ Equipments & WiRe™ App & Web Service.

DAIHAN Smart-Lab™ Equipments

The Smart-Lab™ System is an innovation of laboratory equipments. DAIHAN Smart-Lab™ equipments include two types controller, 4-inch or 7-inch, full-touch screen controller, integrating the most extensively applied Android System and Internet Connectivity. At anytime, anywhere you are, and anyhow, you can check and control all the equipments in your lab using your smart phone, tablet on the web.

No Button. No Knob. Simply Touch!!
DAIHAN SMART line of Lab Equipments with
Full Touch Screen Controller and Internet Connectivity!!

» Features

- 4" or 7" Full Touch Screen TFT LCD
- Ergonomic Design for Intuitive User Interface
- Program Control : Multiple Pattern with Multiple Segments
- Data Recording : On-Screen Temperature Chart
- USB storage for Recorded Data
- Internet Connectivity with WiFi
- Remote Monitoring and Control by WiRe™ Service
- **Storage Management function :**
Record the location of samples in the chamber and find easily afterwards, on screen (for ULT Freezers)
- **Self-Diagnostic Function :**
Automatic Check & Alarm for the Failures in Key Parts
- Real-time Notification of any Failure in the Equipment to the User(s), to the Service Man and to Supplier, through WiRe™ Service



DAIHAN Smart-Lab™ Controller

The DAIHAN Smart-Lab™ Controller is embedded in all of Wisd Smart-Lab™ equipments, which has a lot of useful functions like data recording, Storage management, and self-diagnostic function and so on.

These functions are unique and innovative in the laboratory equipments. And the controller is connected with researcher's Smart phone, tablet or laptop computer through internet, by WiFi, so scientists are able to check all the equipments in the laboratory at any time, and at any place.



» Ergonomic Screen Design

- 4" or 7" Full Touch Screen TFT LCD
- Intuitive User Interface
- No-extra Button, Knob or Switch
- All the main Functions and Parameters are displayed on the Main Screen.
- Simple Touch to set the Temperature, Time, and other parameters



» Automatic Data Recording

- High Resolution Graphic Display
- On-Screen Temperature Chart : Check the Temperature Record Graph on the Screen and could see the Details by Zoom-in & Zoom-out enabled at your finger tips.
- The Recorded Data could be transferred to PC by USB Memory.

No.	Start Segment	End Segment	Count	Delete
01	No Segment !	No Segment !	No Count !	Delete
02	No Segment !	No Segment !	No Count !	Delete
03	No Segment !	No Segment !	No Count !	Delete
04	No Segment !	No Segment !	No Count !	Delete

» Program Control

- Variable Program Setting for Multiple Patterns with Multiple Segments
- Preview window offered to check simply the input patterns
- Very easy to input and edit patterns



» Self Diagnostic Function

- Automatic Check & Alarm for the Failures in Key Parts of the equipment
- Real-time Notification of any Failure in the Equipment to the User(s), to the Service Man and to Supplier, through WiRe™ Service.
- Implemented by LIN and empowered by WiRe™ Service




» Internet Connectivity with WiFi

- Remote Control is available for all the Equipments in the laboratory using Smart Phone, tablet or on the Web.
- Available WiRe™ App & Web Service

Maximum Convenience by Advanced Digital Technology

Smart-Lab™ Equipments



WiRe™ App & Web Service




WiRe™ App & Web service enables gathering, analyzing, and sharing of experimental data, and direct control of the unit is possible as well, anywhere and at any time.

When a defect of the unit occurs, it is processed by self-diagnostic system automatically. If you want to get safe and reliable data comfortably, it's a Smart choice to be Supported by WiRe™ service.

*At anytime, Anywhere you are, and Anyhow, you can check and Control All the Equipment in your Lab **Using Your Smart Phone or on the Web.***

» Features

- Developed by Advanced Web Technology
- Available in all the PCs, Laptops and MACs
- Almost same UI as the Touch Screen Controller on the Equipment is available
- Support all different kinds of Web browser based on HTML Standards
- Offering Android App and Apple iOS App



» WiRe™ Remote Control System

- Anytime, anywhere you are, and anyhow, you can check and control the Freezer in your Lab using your smart phone or web.
- It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet and web without restraint of time and place.



» Device Settings and User Registration

- A single unit can be remote controlled by multiple users
- A single user can remote control multiple units in the Lab
- Smart-Lab™ equipments offer registration of multiple users with simple account information



» Self Diagnostic Function

The Self-Diagnosis System senses a certain defect within the unit, and reports it to the service center through the servers operated by DAIHAN. The Service Center verifies the defects of the unit received, and could offer prompt service.



Display Error Status Reporting

» Smart-Lab™ Equipment lines

ThermoStable™ SIG- / SIF-Series of General Purpose Incubators	P.112~113, 116~117
ThermoStable™ SIR-Series of Low Temperature Incubators	P.126~127
ThermoStable™ SON- / SOF-Series of General Purpose Ovens	P.148~149, 152~153
ThermoStable™ SOV-Series of Vacuum Ovens.....	P.160
ThermoStable™ GC-Series of Growth Chambers	P.48~49, 130~131
SimpleFreez™ & DuoFreez™-Series of Ultra Low Temperature Freezers	P.166~181
SMART Refrigerated Bath Circulators	Coming soon
SMART Shaking Incubators	Coming soon
SMART Temp. & Humidity Chambers	Coming soon
SMART Cold Lab Chambers	Coming soon

DAIHAN ThermoStable™ SOV SMART Vacuum Ovens, Highly Safety Viewing Window, 18.6-/30-/69-Lit., up to 200°C with Smart-Lab™ Controller, 4" Full Touch Screen, Fuzzy-PID Control, Superior Temp. & Vacuum Accuracy, 10~750 mmHg NEW



SMART Vacuum Oven "ThermoStable SOV-30"

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Suitable for Drying, Baking, Conditioning, Curing, Out-Gassing Solids and Liquids, Vacuum Embedding, Moisture Testing, Plating and Aging Test
- Very Easy and Convenient for Vacuuming and Venting, Easy to read Vacuum Gauge
- Highly Safety Viewing Window (Tempered Safety Glass+ Polycarbonate Window)
- Smart-Lab™ Controller
- PL(Product Liability) Insurance
- UL/CUL Certified
- CE Certified

- * RS232C Interface for Remote Monitoring and Controlling with PC
- * Superior Uniformity and Stability
- * High Quality Door Sealing : High Temperature Grade Silicone Molded Gasket
- * Over Temperature & Over Current Protection and Sensor Error Detection
- * Ergonomic Design (Smart-Lab™ Controller)
- * WiRe™ App & Web Service.
- * Automatic Data Recording
- * Self-Diagnostic Function
- * Internet Connectivity with WiFi
- * Variable Program Setting for Multiple Patterns with Multiple Segments
- * The Recorded Data Could be transferred to PC by USB memory
- * Controller SW can be upgraded through USB port.

Specification

Model	ThermoStable SOV-20	ThermoStable SOV-30	ThermoStable SOV-70	
Capacity	18.6 Lit.	30 Lit.	70 Lit.	
Vacuum Range	10 ~ 750 mmHg			
Temp.	Range / Sensor Ambient +5°C to 200°C / PT100			
	Accu./ Uniformity ± 0.5°C at 100°C / ± 2.0°C at 100°C			
Heating Power/Consumption	300W×2ea	300W×2ea	400W×4ea	
Heat-up Time	50 min. to 100 °C, 90 min. to 150 °C			
Timer & Alarm	99hr 59min (with delay / Continuous run.), Error status & Timer-end			
Display	4" Full Touch Screen TFT LCD			
Controller	Smart-Lab™ Controller			
Resolution	Control : ± 1.0°C, Display : ± 1.0°C			
Shelves	Included 2ea Aluminum Shelves			
Safety Circuit	Over Temp. and Over Current Protector, Sensor Error Detector, Leakage Breaker			
Material	Internal : Stainless steel (# 304), External : Powder Coated steel, Window : Tempered Safety Glass 12T, Poly Carbonate 5T for Safety			
Others	Program Function, WiRe™ App & Web Service.			
Dimension (w×d×h)	Interior (mm)	265×290×265	300×330×300	400×435×400
	Exterior (mm)	472×443×542	507×483×577	570×588×677
Vacuum Inlet Valve Size	od. 10 mm			
Net Weight	48 kg	56 kg	82 kg	
Packing Size (w×d×h) (mm) & Gross Weight		572×660×840	710×762×970	780×826×1,030
		67 kg	76 kg	114 kg
Power Consumption	600 W	600 W	1.7 kW	
Power Supply & Cord Plug*	1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			

* Other Specifications are available upon Customer's Request.

Ordering Information

- DH.SWOV03020 Vacuum Oven, 230V, 18.6 Lit., "ThermoStable SOV-20"
- DH.SWOV03030 , 30 Lit., "ThermoStable SOV-30"
- DH.SWOV03070 , 70 Lit., "ThermoStable SOV-70"
- DH.SWOV04020 Vacuum Oven, 120V, 18.6 Lit., "ThermoStable SOV-20"
- DH.SWOV04030 , 30 Lit., "ThermoStable SOV-30"
- DH.SWOV04070 , 70 Lit., "ThermoStable SOV-70"

Spare Aluminum Shelves

- DH.WOV11020 (1) Standard Shelf, Aluminum, "OVS020", w295×d285×h98mm, for 18.6 Lit.
- DH.WOV11030 , "OVS030", w295×d285×h98mm, for 30 Lit.
- DH.WOV11070 , "OVS070", w395×d395×h130mm, for 70 Lit.

Validation (IQ, OQ) Service

- DH.SWOV31020 Validation Service(IQ, OQ), "SOVV020", for "ThermoStable SOV-20"
- DH.SWOV31030 , "SOVV030", for "ThermoStable SOV-30"
- DH.SWOV31070 , "SOVV070", for "ThermoStable SOV-70"



Display of Smart-Lab™ Controller

※ See page 6 ~ 11 for the Feature of Smart-Lab™ System



(1) Standard Shelf, Aluminum (Included)

Cat. No	Description	Articles
---------	-------------	----------

Vacuumer® Precision Vacuum Pumps, “VOP”-Series, 50 ~ 120Lit. with Double-Stage, Direct Drive & Oil Sealed Rotary Type, Gas Ballast Valve Installed

■ Gas Ballast Valve Installed

- * Double Stage Vacuum Pump
- * Over Current Protection
- * Physical and Chemical Apparatus Industries
- * Industrial Machines
- * Vacuum Pumping System

■ Oil Back Flow Prevention Mechanism Installed

- * Easy-to-use Large Finger Grips
- * Installed Filter Prevent Sludge Contamination
- * Vacuum Evaporation Equipment
- * Vacuum Dryer
- * Vacuum Freeze Dryer. Etc

Specification

Model	“VOP-60”	“VOP-100”
Pumping Speed (50/60Hz)	50 / 60 Lit./min	100 / 120 Lit./min
Ultimate Pressure (Pa)	6.7×10^{-2} (6.7×10^{-1})	
Motor Power	Single Phase 220V, 200W, 4poles	Single Phase 220V, 400W, 4poles
Motor Speed(50/60Hz)	1430 / 1730 RPM	
Recommended Oil	COSMA-100	
Oil Capacity(mL)	800	900
Suction Pipe(o.d×i.d)	NW25(KF25)	
Temp. Range	7~40℃	
Weight	17kg	25kg
Dimensions(w×d×h)(mm)	165×411×234	171×464×238

* Per Customer's request various specifications are available.

* Using the Macleod gauge the ultimate pressure is measured. The measurements from the pirani gauge are show ln()

Ordering Information

- V3.VOP060 (1) Vacuum Pump, “VOP-60”, 50 / 60 Lit./min, (50 / 60Hz), 220V
 V3.VOP100 (2) , “VOP-100”, 100 / 120 Lit./min, (50 / 60Hz), 220V

Optional Accessories for VOP-series

- V3.OMT150 (3) Oil Mist Trap for “VOP-60” / “VOP100”
 V3.COS101 (4) Oil for Vacuum Pump VOP-series, 1Lit.
 V3.COS104 (5) , 4Lit.

※ Selection Guide

We Recommend “VOP-60” Vacuum Pump for 18.6- / 30-Lit. Vacuum Oven(P.159~160)
 We Recommend “VOP-100” Vacuum Pump for 70Lit. Vacuum Oven(P.159~160)



(1) VOP-60



(2) VOP-100



(3)



(4) 1Lit.

(5) 4Lit.