Laboratory Solutions



Laboratory Solutions

-	
pH Buffers	
Conductivity S	Standards
Solutions for	ORP and DO
Solutions for	ISE
Maintenance	Solutions



Solutions for Calibration & Care

Comprehensive Range of Indispensable Utilities



Genie in a Bottle a Full Pack of Competence

The determination of pH, conductivity, ion concentration, redox potential and dissolved oxygen are common analyses in most laboratories. The measurement accuracy highly depends on the quality and operation of the solutions used for sensor calibration and maintenance. METTLER TOLEDO takes pride in a long tradition of producing complete measurement systems, including a comprehensive range of top quality solutions.

Top sellers for the most common use



All METTLER TOLEDO calibration and cleaning solutions are available in handy bottles with no-fuss labels. If larger quantities are required, then our 6-packs are the right choice. Electrolytes come with special caps that make refilling your electrode a child's play.

Specialties for specific requirements



Using our sachets guarantees fresh solution for every calibration as well as maximum ease of use. Every sachet box comes with a printed certificate. Rely on METTLER TOLEDO's expertise to meet your specific needs with our expanded offering.

Maximum traceability and compliance



To guarantee maximum traceability, an individual test certificate exists for every calibration solution. Furthermore, compliance with regulations is key for us. All SDS and lables contain information according to GHS (Globally Harmonized System) in local languages. Simply download all required documents from our online database: www.mt.com/buffer





Good Electrochemistry Practice[™] for calibration solutions

The quality of a solution is only guaranteed for unopend bottles that have not yet expired. The following tips shall help optimize the use of calibration solutions after opening and while in use in the laboratory. Fresh calibration solutions minimize the measurement uncertainty and thus optimize the reproducibility of results:

- Note the date on the bottle when the calibration solution is opened for the first time.
- Keep the bottles sealed tightly at all times and use any dispensed calibration solution promptly.
- Never pour dispensed calibration solution back into the original bottle.
- Ensure that bottles containing calibration solution are not contaminated.
- Store calibration solutions at ambient temperatures and avoid direct sunlight.
- Clean your sensors properly prior to calibration and do not calibrate directly in the original bottle (unless there is an instruction to do so).
- Replace calibration solutions that have reached the expiry date or may be contaminated.

www.mt.com/GEP

pH Measurement Precision Starts with a Precise Calibration

pH measurements are only as accurate as the buffer solutions used for calibration purposes. METTLER TOLEDO offers a selection of quality pH buffers to match your specific requirements. No matter if you look for traceable technical buffers or buffers that are certified by an accredited body, you will find the right solution. Maximum precision is guaranteed with our NIST/DIN buffers!



The internationally recognized pH scale is based on standard reference materials (SRM) selected by NIST (National Institute of Standards and Technology, USA). METTLER TOLEDO buffer solutions are traceable to these primary standards. The quality inspection certificate, that is available for every bottle and sachet, guarantees stated values and traceability.

Traceable to an accredited body: DKD certified buffers



The DKD (Deutscher Kalibrierdienst, German Calibration Service) certified buffers are the perfect product for regulated industries. DKD is accepted by other accreditation bodies within the scope of the European Cooperation for Accreditation (EA) and the International Laboratory Accreditation Cooperation (ILAC).

Top Precision: NIST/DIN 19266 buffers



The NIST/DIN pH buffers are specified to three decimal places, e.g. 9.180, offering the highest possible pH precision. A detailed test certificate is available for every bottle and guarantees the tested values and traceability. Calibration cannot get better than this!



Automatically compensated temperature



The pH value of a buffer solution changes with temperature. Several international standards state the pH values at different temperatures in steps of 5 °C. Thus, all METTLER TOLEDO buffer solutions come with a temperature table printed on the label. The same tables are stored in all METTLER TOLEDO pH meters which allows for easy automatic temperature compensation.



Good Electrochemistry Practice[™] for reliable pH calibrations

pH electrodes must be calibrated regularly; at least once a day prior to measurements and after every cleaning, regeneration or long-term storage. The calibration adjusts slope and offset of an electrode to their true values. Since an electrode is characterized by both its slope and its offset, it is advisable to do a minimum of a two point calibration for reliable measurements and better precision. For a wider measurement range three or more calibration points are recommended. Important is that the following measurements are within the calibration range.

Optimal Care For Hard-Working pH Electrodes

The pH electrode is the main actor of the measurement, requiring, therefore, special care and attention. METTLER TOLEDO provides the necessary maintenance solutions. For an easy start the All-in-One kits offer an optimal selection of solutions for calibration and maintenance. Your electrode stays ready for whatever task is due. Reliable pH results guaranteed!

The right electrolyte for every application



Thanks to the ARGENTHAL[™] reference system most METTLER TOLEDO pH sensors are filled with 3 mol/L KCI electrolyte solution. There is no more risk of contaminating your samples with silver ions. For measurements in nonaqueous samples or samples with low ionic strength a special bridge electrolyte is needed. Always make sure that your electrode is completely filled with clean electrolyte solution.

Maintenance solutions for peace of mind



When rinsing with deionized water is not sufficient, a special cleaning solution can be used to remove sample residues. Depending on the contamination, Pepsin-HCl or Thiourea is recommended. The InLab[®] storage solution provides optimal conditions for sensors during the time between measurements, be it short- or long-term storage. Find more helpful maintenance tips on **www.electrodes.net.**

Quick and easy performance verification



The VPac[™] pH verification kit provides a quick and easy performance verification of your pH system. Measure two test solutions of unknown pH value and get an immediate online evaluation of your results including certificate! Never has it been easier to verify your system after installation or any setup change or simply for extra peace of mind.

www.mt.com/pH-VPac



pH Buffers	pH value at 25 °C	Order number 250 mL	Order number 6 x 250 mL	Order number 30 sachets 20 mL
Technical pH	2.00	51350002	51350016	30111134
buffer solutions	4.01	51350004	51350018	51302069
	7.00	51350006	51350020	51302047
	9.21	51350008	51350022	51302070
	10.00	51350010	51350024	51302079
	11.00	51350012	51350026	30111135
	Rainbow bottles I (3 x 2 bottles 250 mL 4.01 / 7.00 / 9.21)		30095312	
	Rainbow bottles II (3 x 2 bottles 250 mL 4.01 / 7.00 / 10.00)		30095313	
	Rainbow sachets I (3 x 10 sachets 20 mL 4.01 / 7.00 / 9.21)			51302068
	Rainbow sachets II (3 x 10 sachets 20 mL 4.01 / 7.00 / 10.01)			51302080
NIST/DIN pH buffer solutions	4.006	51350052		30111136
	6.865	51350054		30111137
	9.180	51350056		30111138
	10.012	51350058		30111139
Certified pH	4.01	51350032	51350042	
buffer solutions	7.00	51350034	51350044	
	9.21	51350036	51350046	
	10.00	51350038	51350048	

Electrolytes for reference electrodes	Order number 25 mL	Order number 250 mL	Order number 6 x 250 mL	Order number 6 x 30 mL
KCI-solution 3 mol/L for ARGENTHAL [™] reference systems	51343180	51350072	51350080	
KCI-solution 3 mol/L, AgCI saturated, for Ag/AgCI reference systems	51343184	51350074	51350082	
FRISCOLYT-B [®] , for measurement at low temperature and for media with organic compounds (oil, proteins etc.)	51343185	51350076	51350084	
LiCl Solution 1 mol/L in ethanol, for measurement in non-aqueous media				51350088

Maintenance solutions	Order number 250 mL	Order number 6 x 250 mL	Order number 25 mL
Pepsin-HCI for cleaning junctions with protein contamination. Treatment time about 1 h.	51350100	30045061	
Thiourea solution for cleaning junctions with silver sulfide contamination. Treatment until discoloration.	51350102	30045062	
Reactivation solution for regeneration of glass electrodes. Treatment time about 1 min.			51350104
InLab [®] storage solution for pH and ORP electrodes	30111142		
pH All-in-One Kit I (pH buffer 4.01 / 7.00 / 9.21, 3 mol/L KCl, cleaning solution, storage solution)		30095314	
pH All-in-One Kit II (pH buffer 4.01 / 7.00 / 10.00, 3 mol/L KCI, cleaning solution, storage solution)		30095315	
	Order number		

	2 x 100 mL
VPac [™] pH verification kit	30090849

Conductivity Standards and More It's all About Correct Handling

Depending on the type of sensor, conductivity standards are used for calibration or verification. Low conductivity standards require special handling and are mostly used for verification purposes. In this case the cell constant has been determined with a standard of higher conductivity or is stated on the sensor certificate. Conductivity standards are available in sachets to guarantee fresh solution for every calibration as well as maximum ease of use.



Low conductivity standards – the influence of air

Conductivity standards are directly affected by the influence of carbon dioxide (CO_2) when in contact with air. Thus, especially the lower conductivity standards have a limited lifespan. Measurements of samples with conductivity lower than 10 µS/cm need a special procedure, such as protection with inert gas or the usage of a flow cell. These standards are only intended for verification and not for calibration.



Temperature dependence

A slight change in temperature usually has a big impact on the conductivity value of a standard solution. A table on every bottle label indicates the conductivity values at the most common measurement temperatures. During calibration the meter automatically references to this table for temperature compensation. If possible, calibration and measurements should be performed at the same temperature.

Order number 250 mL	Order number 6 x 250 mL	Order number 30 sachets 20 mL
30090847		
30094617		
51300169		
51302153		
51300170		
51350092	51350096	51302049
51350094	51350098	51302050
	250 mL 30090847 30094617 51300169 51302153 51300170 51350092	250 mL 6 x 250 mL 30090847 30094617 51300169 51302153 51300170 51350096

* Maximum storage: 1 month

**Maximum storage: 3 months



Redox Buffer Solutions and Tablets for Dissolved Oxygen Sensors



Redox buffer solutions for verification purposes

Redox buffer solutions are used for verification of all common redox sensors. They are not used for calibration purposes. Similar to other solutions they are temperature dependent. It is therefore important to know the measurement temperature of the buffer. A table on every bottle label indicates redox values at different temperatures.



Zero oxygen tablets

Zero oxygen tablets make the preparation of a solution with zero oxygen content very easy. This solution can be used for calibration, verification or conditioning purposes if measurements are performed at low dissolved oxygen levels.

Redox buffer solutions	E (Ag/AgCI) 25 °C	Order number 250 mL	Order number 6 x 250 mL	Order number 6 x 30 mL
	220 mV, pH 7 (U _H = 427 mV)	51350060	51350062	
	468 mV, pH 0.1 (U _H = 675 mV)			51350064
DO Accessories		Order number		
Zero oxygen tablets (20 pcs.)		51300140		

Solutions for Ion-Selective Electrodes The Proper Mix for Accurate Results

Measuring with ion-selective electrodes (ISE) is the easiest and most affordable way to determine ion concentration. However, ISEs require careful handling and the use of the correct solutions. METTLER TOLEDO offers all solutions that are needed for successful ion measurements.

Ready to use ion calibration standards



High precision ion calibration standards can be ordered at concentrations of 1000, 100 and 10 mg/L (ppm). In case a lower concentration is needed it can easily be prepared with serial dilution that is explained in the electrode manual.

Ionic Strenght Adjustors for high repeatability





In all analytical procedures using an ISE, the correct amount of ISA (lonic Strenght Adjustor) must be added to all samples and standards prior to measurement or calibration. This solution ensures that samples and standards have similar and constant ionic strength. Instructions for type and amount of ISA can be found in the manual of every ISE.



It is important to fill the reference electrode of every ISE with the recommended electrolyte solution. The correct electrolyte will minimize junction potentials and provide optimum temperature and time response. The electrolyte must be refilled or replaced regularly in order to achieve good electrode performance. For more information please refer to the electrode manual.



Solutions for perfectION $^{\scriptscriptstyle \rm M}$ combined ISE

Reference electrolyte solutions	Order number 5 x 60 mL
Ion Electrolyte A (calcium, fluoride, sulfide)	51344750
Ion Electrolyte B (chloride, cyanide, lead, silver/sulfide)	51344751
Ion Electrolyte C (silver)	51344752
Ion Electrolyte D (copper, iodide)	51344753
Ion Electrolyte E (potassium)	51344754
Ion Electrolyte F (nitrate)	51344755

ISA solutions	Order number 475 mL	Order number 3790 mL
ISA solid state ISE (chloride, copper, iodide, silver)	51344760	
Calcium ISA	51344761	
Potassium ISA	51344762	
Nitrate ISA	51344763	
Nitrate ISS (for suppressing interference)	51344764	
Fluoride TISAB II with CDTA		51344765
Fluoride TISAB III with CDTA (concentrate)	51344766	

Solutions for DX series ISE half-cells

Bridge electrolyte	Order number 25 mL	Order number 250 mL	Order number 6 x 250 mL
1 mol/L KNO ₃	51343182	51350078	51350086
3 mol/L KCl	51343180	51350072	51350080
1 mol/L KCl	51343181		

ISA solutions

TISAB 3, for fluoride determinations	51350106
0.9 mol/L Al ₂ (SO ₄) ₃	51350108

Ion calibration standards

	Order number 500 mL	Order number 500 mL	Order number 500 mL
	1000 mg/L	100 mg/L	10 mg/L
Silver ISE standard solution	51344770		
Calcium ISE standard solution	51344771	30090855	30090856
Chloride ISE standard solution	51344772	30090853	30090854
Cyanide ISE standard solution	51344773		
Copper ISE standard solution	51344774		
Fluoride ISE standard solution	51344775	30090851	30090852
lodide ISE standard solution	51344776		
Potassium ISE standard solution	51344777		
Sodium ISE standard solution	51344778	30090857	30090858
Ammonium ISE standard solution	30090859	30090860	
Nitrate ISE standard solution	51344779		
Lead ISE standard solution	51344780		
Sulfide ISE standard solution	51344781		

Discover the Safest Path to the Top with Good Electrochemistry Practice[™]

Various factors can affect your pH, redox, conductivity, dissolved oxygen or ion measurements. Take 5 minutes to localize your risks and get the neccessary support:

statication staticatii statication statication statication statication statica

www.mt.com/BuffersAndMore

www.mt.com/GEP

For more Information

Mettler-Toledo AG, Analytical CH-8603 Schwerzenbach, Switzerland Tel. +41 22 567 53 22 Fax +41 22 567 53 23

Subject to technical changes © 03/2014 Mettler-Toledo AG, 30126698 Marketing pH Lab / MarCom Analytical ISO 9001 regissured

Environmental management system according to ISO 14001.

and testing according to ISO 9001.

Quality certificate. Development, production



"European conformity". The CE conformity mark provides you with the assurance that our products comply with the EU directives.