DAIHAN-brand[®] Digital Precise Water Bath, "WB", 6-/11-/22-Lit., up to 100°C, ±0.1°C, CE-Certi. with Digital Fuzzy Control System, 1×Stainless steel Flat Lid, and Back Light LCD, with Certi. & Traceability 정밀 항온수조, 디지털 퍼지 제어 시스템, Stainless steel Flat Lid 포함



Digital Precise Water Bath, 11Lit. "WB-11"



(1) Flat Lid, Stainless steel (Included)



(3) "BFC106" Lid with Concentric Ring Cover Sets for 6Lit.



(5) "BFC122" Lid with Concentric Ring Cover Sets for 22Lit. Bath



(2) Dome Lid, Stainless steel



(4) "BFC111" Lid with Concentric Ring Cover Sets for 11Lit.



(6) Test Tube Rack, Stainless steel

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
 Patented Jog-Shuttle Control System and Digital LCD

Description

- with Back-Light Function PL(Product Liability) Insurance
- * User's Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature (\pm 10.0°C).
- * 304 Stainless steel Bath & Lid for Superior Durability & High Thermal Efficiency/Constant.
- * Digital Timer Function(99hr 59min) : 'Delay Time' and 'Operation Time' can be Controlled.
- * Storage Function for the Set Values of Temperature and Timer.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled).

Specification

Capacity & Models		6 Lit. "WB-6"	11 Lit. "WB - 11"	22 Lit. "WB - 22"		
Temp. Range & Accuracy		Ambient+5℃~ 100℃, ±0.1℃				
Temp. Resolution		0.1 °C − Display, 0.1 °C − Control				
Temp. Uniformity & Probe		±1.0°C, PT100				
Heating Power		800 W	1 kW	2 kW		
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end				
Display		Digital LCD with Back-Light Function				
Controller		Digital Fuzzy Controller with Jog - Shuttle Switch (Turn + Push)				
Safety Circuit		Over-Temp. / -Current protection, Sensor Error Detection				
Material		Stainless steel (#304)-Internal, Powder Coated steel-External				
Lid		included Stainless steel Flat Lid, Optional Dome Lid				
Drain Valve		No	Mounted(ϕ 9.5~11mm)			
Dimension (w×d×h)	Effective Space	300 imes 155 imes h110mm	$302 \times 240 \times h110$ mm	500 imes 295 imes h110mm		
	Internal	$300\!\times\!155\!\times\!h150mm$	$302 \times 240 \times h150mm$	500 imes 295 imes h150mm		
	External	$325 \times 180 \times h260mm$	$325 \times 265 \times h260$ mm	535×330×h260mm		
Packing Size (mm) and Gross Weight (kg)		350×210×h310mm, 8 kg	350×300×h310mm, 10 kg	570×360×h310mm, 14 kg		
Power Supply*		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz				

Ea Price ₩

* Other Specifications are available upon Customer's Request.

Orc	lering	Int	form	nation
Cit	ici nig			auon

DH.WB000106 DH.WB000111 DH.WB000122	Water Bath, Digital, 230V, 6 Lit., "WB-6" , 11 Lit., "WB-11" , 22 Lit., "WB-22"	427,600 456,440 527,780			
DH.WB000206 DH.WB000211 DH.WB000222	Water Bath, Digital, 120V, 6 Lit., "WB-6" , 11 Lit., "WB-11" , 22 Lit., "WB-22"	427,600 456,440 527,780			
Spare Flat Lids, Stainless steel					
DH.WB010106 DH.WB010111 DH.WB010122	(1) Flat Lid, for 6 Lit., w325×d176mm, " BF106 " , for 11 Lit., w325×d265mm, " BF111 " , for 22 Lit., w530×d325mm, " BF122 "	20,840 22,640 29,940			
Optional Lids and Tube Rack, Stainless steel					
DH.WB010206 DH.WB010211 DH.WB010222	(2) Dome Lid, for 6 Lit., 324 × 176 × h100mm, "BD206" , for 11 Lit., 324 × 263 × h100mm, "BD211" , for 22 Lit., 533 × 328 × h100mm, "BD222"	132,430 146,430 166,830			
DH.WB030106	(3) Lid with 2×Concentric Ring Cover Sets, for 6 Lit., Stainless steel, 330×182×h20mm, "BFC106"	64,150			
DH.WB030111	(4) Lid with 2 × Concentric Ring Cover Sets, for 11 Lit. , Stainless steel, 329 × 269 × h20mm, " BFC111 "	78,470			
DH.WB030122	(5) Lid with 6 × Concentric Ring Cover Sets, for 22 Lit. , Stainless steel, 539 × 334 × h20mm, " BFC122 "	137,890			
DH.WB020106	(6) Tube Rack, for 6 Lit ., 36Hole $\times \phi$ 17mm, " BTH106 ", Stainless steel, 220 \times 125 \times h78.5mm	35,830			
DH.WB020206	(6) Tube Rack, for 6 Lit ., 15Hole × <i>Φ</i> 32mm, " BTH206 ", Stainless steel, 220 × 125 × h78.5mm	35,830			
DH.WB020111	(6) Tube Rack, for 11 Lit ., 63Hole × <i>Φ</i> 17mm, "BTH111" , Stainless steel, 220 × 205 × h78.5mm	42,780			
DH.WB020211	 (6) Tube Rack, for 11 Lit., 25Hole × <i>Φ</i> 32mm, "BTH211", Stainless steel, 220 × 205 × h78.5mm 	42,780			
	Validation (IQ, OQ) Service	Ea. Price, ₩			
DH.WB040106	Validation Service(IQ, OQ), "BV106", for "WB-6"	문의요함			
DH.WB040111 DH.WB040122	, "BV111", for "WB-11" , "BV122", for "WB-22"	문의요함 문의요함			