

**DAIHAN-brand® High-performance Dry Bath, Heating, Dual-type “HB-96D”, up to 150°C, ±0.1°C**  
with Digital Fuzzy Control, Molded Heater, Modular Blocks, Acrylic Lid, Touch-button Controller, with Certi. & Traceability



Dry Bath, Heating, Dual Control System  
“HB-96D-set”

- with Certificate & Traceability :  
Controlled by Serial Number, Certificate, Delivery-information,  
and Traceable Data Base System
- Designed for Multi-functional Dry Block Heater
- Dual Controller & Block
- **Patented Jog-Shuttle Control System** and Digital LCD  
with Back-Light Function
- CE Certified                      ■ **PL (Product Liability) Insurance**

- \* Ideal for Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures
- \* Precise Temperature Control : ±0.1°C at 37°C, ±0.2°C at 90°C, ±0.4°C at 135°C
- \* User Self-Compensation Function to offset between Actual and Displayed Temperature. (±10.0°C)
- \* Over Temperature & Over Current Protection and Sensor Error Detection
- \* Digital Timer Included : up to 99hr 59min (Operation Time + Delay Time)
- \* Storage Function of Temperature and Timer setting values
- \* Alarm Function : Error Status and Timer-end
- \* Various kind of Blocks are offered (Optional),                      \* Modular Block
- \* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- \* The Block Extraction Tongs and Included for Easy Handling

### Specification

Models	“HB-96D”
Controller	Dual Digital Fuzzy Controller
Capacity (1.5/2ml Tube)	2×(48 Holes×1.5/2ml Tube)
Heater	300W (Molded Heater)
Dimension	Block (mm)
	Unit (mm)
Temp. Range, & Accuracy	Ambient +5°C to 150°C, ±0.1°C at 37°C, ±0.2°C at 90°C, ±0.4°C at 135°C
Heat-Up Speed & Sensor	Approx. 5°C/min, PT100
Safety Circuit	Over Temp. and Over Current Protector, Sensor Error Detector
Timer with Alarm	99hr 59min, Alarm (Error Status and Timer - end)
Display	Digital LCD with Back-Light Function
Material	Internal & External
	Block
Included Standard Acce.	2×Block “BLC548” for 1.5/2ml Tubes, and Block Tongs “BLE101”
Packing Size	400×270×h200 mm
Gross Weight	6 kg
Power Consumption	600 W
Power Supply & Cord Plug *	1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz

\* Other Specifications are available upon Customer's Request.

### Ordering Information

#### HB-96D with 1.5/2ml Block-set

DH.WHB00199 Dry Bath, Heating, 230V, “HB-96D-set”, with 1.5 / 2ml Block  
DH.WHB02199 , 120V, “HB-96D-set”, with 1.5 / 2ml Block

#### HB-96D only unit

DH.WHB00198 Dry Bath, Heating, 230V, “HB-96D-unit”, without Block  
DH.WHB02198 , 120V, “HB-96D-unit”, without Block

### Included Block in the Package-Set



(3) “BLC548”  
2×1.5/2ml Centrifuge Tube×48 holes

<Continued on next page>

Cat. No	Description	Articles
---------	-------------	----------

<Continued from... DAIHAN-brand® High-performance Dry Bath, Heating, Dual-type "HB-96D">

Ordering Information

Optional Interchangeable Blocks

- DH.WHB510596 (1) Block, "BLC596", for 0.2ml PCR Tube, 96 holes
- DH.WHB510570 (2) Block, "BLC570", for 0.5ml Centrifuge Tube × 63 holes
- DH.WHB511548 (3) , "BLC548", for 1.5/2ml Centrifuge Tube × 48 holes
- DH.WHB521248 (4) Block, "BLT248", for φ 12mm Test Tube, 48 holes
- DH.WHB521335 (5) , "BLT335", for φ 13mm Test Tube, 30 holes
- DH.WHB521535 (6) , "BLT535", for φ 15mm Test Tube, 30 holes
- DH.WHB521628 (7) , "BLT628", for φ 16.5mm Test Tube, 24 holes
- DH.WHB521828 (8) , "BLT828", for φ 18mm Test Tube, 24 holes
- DH.WHB531515 (9) Block, "BLC515", for 15ml Conical Tube, 15 holes
- DH.WHB535008 (10) , "BLC008", for 50ml Conical Tube, 8 holes

DH.WHB500101 Spare Block Extraction Tongs, "BLE101"



Validation (IQ, OQ) Service

DH.WHB30198 Validation Service(IQ, OQ), "DBV198", for "HB-96D"



Dry Bath, Heating, Dual Control System "HB-96D-set"

The Details of 10 Optional Blocks



(1) "BLC596"  
0.2ml PCR Tube × 96 holes



(2) "BLC570"  
0.5ml Centrifuge Tube × 63 holes



(3) "BLC548"  
1.5/2ml Centrifuge Tube × 48 holes  
Included in the Package-Set



(4) "BLT248"  
φ 12mm Tube × 48 holes



(5) "BLT335"  
φ 13mm Tube × 30 holes



(6) "BLT535"  
φ 15mm Tube × 30 holes



(7) "BLT628"  
φ 16.5mm Tube × 24 holes



(8) "BLT828"  
φ 18mm Tube × 24 holes



(9) "BLC515"  
15ml Conical Tube × 15 holes



(10) "BLC008"  
50ml Conical Tube × 8 holes