

Cat. No	Description	Articles
---------	-------------	----------

DAIHAN-brand® Programmable Digital Rockers, "RK", Digital Feedback Control, 5~50 rpm, with Certi. & Traceability with 290×200mm or 300×300mm Platform, Smooth Start & Wave Motion, Tilt Angle $\pm 6^\circ / \pm 9^\circ$

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- High Precision Digital Speed Control : 1rpm
- Patented Jog-Shuttle Control System
- CE Certified ■ PL(Product Liability) Insurance

- * Gel Staining, Destaining, Blotting
- * Cell Culture, Hybridization, Washing
- * High Precision Digital Speed Control : 1rpm
- * Programmable Operation Mode Supported : 10 Steps and 6 Programs
- * Digital Timer : Max. 99hr 59min
- * Digital LCD with Back-Light Function for Timer, RPM and Power Bar-Graph
- * Smooth Start and Change of Rocking Speed
- * Overload Protection Function for Motor
- * Storage Function for Set Values : RPM and Timer
- * Alarm Function : Error Status and Timer-end
- * Stainless steel Platform Comes with Gray Nonslip Rubber Pad
- * Stackable Platform with Pad to Double Mixing Area (Optional)
- * Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

Specification

Model		Platform Size : 290×200 mm "RK-1D"	Platform Size : 300×300 mm "RK-2D"
Dimensions (w×d×h)	Platform (mm)	290×200	300×300
	Overall (mm)	290×240×150	300×330×150
Speed Range (rpm)		5 ~ 50	
Angle (degree)		± 6 or ± 9 (± 9 as Factory Default)	
Motor		35W Geared DC Motor	
Program Function		10 Steps, 6 Programs	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch(Turn+Push)	
Safety Circuit		Overload Protector for motor	
Platform		Included Platform(290×200mm)	Included Platform(300×300mm)
Others		CE Certified, PL(Product Liability) Insurance Smooth Start and Change of the Stirring Speed for Stable Running Storage Function (RPM, Timer) Locking Mode (Jog-Shuttle Input Disabled) Alarm (Error Status and Timer-end)	
Packing Dimension (w×d×h) (mm)		450×350×250	
Gross Weight		7 kg	8 kg
Power Consumption		35 W	
Power Supply & Cord Plug*		Free Voltage (1 Phase, AC 80V~240V, 50/60Hz)	

* Other Specifications are available upon Customer's Request.

Ordering Information

with 290×200 Platform/Rubber Pad, Programmable Digital Rocker

DH.WRK01010 Digital Rocker, 230V, "RK-1D"

DH.WRK02010 , 120V, "RK-1D"

with 300×300 Platform/Rubber Pad, Programmable Digital Rocker

DH.WRK01020 Digital Rocker, 230V, "RK-2D"

DH.WRK02020 , 120V, "RK-2D"

Optional Second Platform

DH.WRK05011 Second Platform, "RK5011", for "RK-1D", w290×d200, with 4 Pillars(L120mm)

DH.WRK05021 , "RK5021", for "RK-2D", w300×d300, with 4 Pillars(L120mm)



Programmable Digital Rocker with Platform "RK-1D"



Programmable Digital Rocker, with Second Platform(Optional), "RK-1D"