


InLab®605

dissolved oxygen



Order number 51340291

former order number -

Specifications

measurement range 0 ... 200 % 0 ... 20 mg/L

compatible with SevenGo pro™ SG6 and
SevenGo Duo pro™ SG68

temperature range 0 ... 60

temperature probe NTC 22 kΩ

shaft material PPS

shaft length 120 mm

shaft diameter 12 mm

connector

electrode head -

cable length 1.8 m

connector BNC/RCA (cinch)

Storage

dry

Special features

IP67

fixed cable

membrane is protected by a steel mesh

shaft made of glass fiber reinforced PPS

Choice of accessories

Membrane kit: 51340293 (containing 3

membrane bodies and 25 mL of electrolyte)

Electrolyte: 51340294