


LE407 - DIN

2-in-1

pH



Order number 51340341
former order number -

Specifications
pH range 0 - 14
zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C
temperature probe -

shaft material POM
membrane glass U
membrane resistance < 250 MΩ

reference system Ag / AgCl
reference diaphragm Nylon
reference electrolyte gel

bridge electrolyte -

electrode head -
cable length 1.0 m
connector DIN

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage 3mol/L KCl

Special features
Unbreakable, chemically resistant POM shaft
Two nylon diaphragms

Scale 1:1