



Order number 51340348

former order number -

Specifications

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C

temperature probe NTC30K

shaft material Glass

membrane glass HA

membrane resistance < 600 MΩ

reference system Ag / AgCl

reference diaphragm Ceramic

reference electrolyte 3mol/L KCl
AgCl saturated

bridge electrolyte -

electrode head LE cap

cable length 1.0 m

connector BNC / Cinch

cell constant (cm^{-1})

cell type -

indicating metal -

Storage 3 mol/L KCl

Special features

Fast measurement in aqueous solution

Integrated temperature probe

Refillable reference electrolyte